



# SPERA™ COVID-19 Ag Test

## Instructions For Use

### For Use with SPERA™ COVID-19 Ag Test

For Use under the Emergency Use Authorization (EUA) only

For *in vitro* diagnostic use only

For prescription use only

For use with kit provided nasal swabs

### Intended Use

The SPERA™ COVID-19 Ag Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigen from SARS-CoV-2 in direct anterior nasal swabs from individuals suspected of COVID-19 by their healthcare provider within the first five (5) days of symptom onset. Testing is limited to laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C §263a, that meet the requirements to perform moderate, high, or waived complexity tests. This test is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

The SPERA™ COVID-19 Ag Test does not differentiate between SARS-CoV or SARS-CoV-2 viruses.

Results are for the identification of SARS-CoV-2 nucleocapsid protein. Antigen is generally detectable in anterior nasal swabs during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Positive results do not rule out bacterial infection or co-infection with other viruses. Additional confirmatory testing with a molecular test for positive results may be necessary if there is a low likelihood of SARS-CoV-2 infection, such as in individuals without known exposures to SARS-CoV-2 or residing in communities with low prevalence of infection. The agent detected may not be the definite cause of disease. Laboratories within the United States and its territories are required to report all results to the appropriate public health authorities.

Negative results should be treated as presumptive and may be confirmed with a molecular assay, if necessary, for patient management. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's respiratory exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19.

The SPERA™ COVID-19 Ag Test is intended for use by medical professionals or operators who are proficient in performing tests in point of care settings. The SPERA™ COVID-19 Ag Test is only for use under the Food and Drug Administration's Emergency Use Authorization.

### Summary and Explanation of the Test

Coronaviruses are a large family of viruses which may cause illness in animals or humans. SARS-CoV-2 is an enveloped, single-stranded RNA virus of the genus. SARS-CoV-2, also known as the COVID-19 virus, was first identified in Wuhan, Hubei Province, China December 2019. The WHO declared that COVID-19 was a pandemic on March 11, 2020, and human infection has spread globally, with hundreds of thousands of confirmed infections and deaths.

The median incubation time is estimated to be 5.1 days with symptoms expected to be present within 12 days of infection. The symptoms of COVID-19 are similar to other viral respiratory diseases and include fever, cough and shortness of breath.

The SPERA™ COVID-19 Ag Test is a rapid lateral flow immunoassay for the qualitative detection of nucleocapsid protein antigen to SARS-CoV-2 in anterior nasal swab specimens directly placed in sample collection tube containing the extraction reagent buffer supplied in the test kit. The SPERA™ COVID-19 Ag Test is supplied with all components necessary to complete a test for SARS-CoV-2.

### Principles of the Test

The SPERA™ COVID-19 Ag Test utilizes immunochromatographic technology that uses antibodies to detect SARS-CoV-2 nucleocapsid protein from anterior nasal swab specimens. Antibodies specific for SARS-CoV-2 as well as control antibodies are immobilized as two distinct lines on the nitrocellulose membrane.

The anterior nasal swab collected from the patient is placed into the sample collection tube to release the specimen containing viral particles from the swab. Three drops of the sample suspension solution are then added to the sample well on the test device. Results are visually interpreted 15 minutes following addition of sample to the device.

### Reagents and Materials Provided

Contents Name	Quantity	Description
Individually packaged SPERA™ COVID-19 Ag Test Devices REF: XH-100-110	10	Anti-SARS-CoV-2 mouse monoclonal antibody immobilized on nitrocellulose membrane, and anti-SARS-CoV-2 mouse monoclonal antibody bonded latex dried in the pad
SPERA™ COVID-19 Ag Test Sample Collection Tubes REF: XH-100-120	0.4 mL per tube x 10 tubes x 5 tubes per bag x 2 bags	Sample collection tube with a buffer solution containing a surfactant and 0.08% (w/v) of sodium azide as a preservative
SPERA™ COVID-19 Ag Test Sample Collection Caps REF: XH-100-130	10	Filter tops for dropping sample suspension buffer containing suspended sample into the device
Sterile sample collection swabs REF: XH-100-140	10	Individually packaged sterile nasal swabs for specimen extraction
Sample tube stand	1	Disposable paperboard tube stand, assembly required before use
Package Insert	1	Instructions for use
Quick Reference Guide	1	Quick reference instructions

### Reagents Materials Required, But Not Provided

Contents Name	Quantity Required	Description
Timer or stopwatch	1	Device to reliably keep track of time
SPERA™ COVID-19 Ag Test Positive Control Swab REF: XH-100-200	1 per new lot of kit or per new operator.	Swab coated with non-infectious recombinant SARS-CoV-2 nucleocapsid protein antigen.  Purchase the SPERA™ COVID-19 Ag Positive Control Swab (REF: XH-100-200) by contacting: support@xtravahealth.com.
SPERA™ COVID-19 Ag Test Negative Control Swab REF: XH-100-205	1 per new lot of kit or per new operator	Sterile blank nasal swab  Purchase the SPERA™ COVID-19 Ag Negative Control Swab (REF: XH-100-205) by contacting: support@xtravahealth.com.

## Warnings and Precautions

- For in vitro diagnostic use.
- For prescription use only.
- For use with kit provided nasal swabs.
- This product has not been FDA cleared or approved; this product is authorized by FDA under an EUA for use by authorized laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C §263a, that meet the requirements to perform moderate, high, or waived complexity tests. This product is intended to be authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.
- The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostic tests for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b) (1), unless the declaration is terminated or authorization is revoked sooner.
- Federal Law restricts this test to sale by or on the order of a licensed practitioner (U.S. only).
- This product has been authorized only for the detection of SARS-CoV-2 antigen, not for any other viruses or pathogens.
- Laboratories within the United States and its territories are required to report all results to the appropriate public health laboratories.
- Do not use specimens stored in media other than SPERA COVID-19 Ag Specimen Buffer.
- Do not store or transport specimens in viral transport media.
- Treat all specimens as potentially infectious. Follow universal precautions when handling samples, this kit, and its contents.
- Use of nitrile or latex gloves is recommended for handling patient specimens.
- Do not use the kit contents beyond the expiration date.
- Do not store product in direct sunlight.
- Do not use product if it has been frozen.
- Leave the test device sealed within its foil pouch until just prior to use.
- Do not use if pouch is damaged or open.
- Do not use a test device that appears damaged or has been dropped after opening
- This product is for single use only. Do not re-use.
- Inadequate or inappropriate sample collection, storage, and transport may yield false test results.
- To obtain accurate results, the Quick Reference Guide must be followed. The SPERA™ COVID-19 Ag Test Instructions for Use is available at [www.xtravahealth.com/spera-product-documentation](http://www.xtravahealth.com/spera-product-documentation).
- Do not use a nasal swab that is not provided with the kit.
- Do not transfer sample from the sample collection tube to the test device without the sample collection cap in place.
- Interpretation of the test before 15 minutes or after 30 minutes could yield false results. Do not interpret the test before 15 minutes or past 30 minutes.
- Reagents contain sodium azide, which is harmful if inhaled, swallowed, or exposed to skin. Contact with acids produces a very toxic gas. If there is contact with skin, wash immediately with plenty of water. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.
- Dispose of used and unused kit contents as biohazardous waste according to federal, state, and local regulatory requirements.
- Wear appropriate personal protective equipment (i.e., clothing, gloves, eye/face protection) when handling the contents of this kit.

## Hazardous Ingredients in the Sample Buffer

Chemical Name/CAS	GHS Code for each Ingredient	Concentration
Sodium Azide / 26628-22-8	Acute Tox. 2(Oral), H300 Acute Tox. 1 (Dermal), H310  MSDS: <a href="https://www.xtravahealth.com/spera-safety-data-sheet">https://www.xtravahealth.com/spera-safety-data-sheet</a>	0.08%

The solution in the tube contains a hazardous ingredient (see table above). If the solution contacts the skin or eye, flush with plenty of water. If irritation persists, seek medical advice <http://www.poisson.org/contact-us> or 1-800-222-1222.

## Kit Storage and Stability

The SPERA™ COVID-19 Ag Test kit should be stored at 2-30°C away from direct sunlight. Do not freeze. The SPERA™ COVID-19 Ag Test kit is stable until the expiration date printed on the outside of the outer packaging, which is currently 12 months from the date of manufacture. If the kit is stored at less than 15°C, ensure all test components are placed at room temperature (15-30°C) at least 30 minutes prior to use. Do not open kit contents until immediately prior to use.

## Quality Control

### Procedural Control

The SPERA™ COVID-19 Ag Test has a built-in procedural control contained within the test device. The blue line at the “control” position of the device will always appear if the sample flows properly and the reagents are working.

### External Positive and Negative Controls

The SPERA™ SARS-CoV-2 Positive Control Swab and Negative Control Swab are available to be provided separately. The external positive and negative controls will validate the entire test. External controls shall be evaluated:

- With every new operator
- Once with each new shipment received (provided that each different lot received in the shipment is tested)
- When problems with testing are suspected or identified
- As necessary to conform with local, state and/or federal regulations, accrediting requirements, or your lab's standard quality control procedures

If the blue control line does not appear or the correct results are not obtained with the external positive and/or negative controls then discard the test and repeat with new components. If the problem persists upon the repeat test, report the problem to Xtrava Health technical support at : 1-888-987-2821.

### Limitations

- This test detects both viable (live) and non-viable SARS-CoV and SARS-CoV-2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- Viral transport media (VTM) should not be used with this test.
- This test is not for use in at-home testing settings.
- False negative results may occur if the swab testing was performed with no swirl mixing or no squeezing in the sample tube.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected or handled improperly.
- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative test results are not intended to rule out other non-SARS viral or bacterial infections.
- Do not use this product in settings with high humidity. False negative results are likely to occur in high humidity conditions.
- Negative results, from patients with COVID-19 symptoms should be treated as presumptive and confirmation with a molecular assay, if necessary, for patient management, may be performed.
- Interpretation of the test after 30 minutes could yield false results. Do not interpret the test past 30 minutes.
- The performance of this test was established based on the evaluation of a limited number of clinical specimens collected between 06/04/2021 and 07/22/2021. The clinical performance has not been established in all circulating variants but is anticipated to be reflective of the prevalent variants in circulation at the time and location of the clinical evaluation. Performance at the time of testing may vary depending on the variants circulating, including newly emerging strains of SARS-CoV-2 and their prevalence, which change over time.
- False positive results may occur, particularly in areas with low numbers of COVID-19 infections and individuals without known exposure to COVID-19 and confirmation with a molecular assay may be considered.
- Positive and negative predictive values are highly dependent on COVID-19 prevalence of disease. False positive test results are more likely during period of low activity when prevalence is moderate or low.

## Conditions for Intended Authorization

The SPERA™ COVID-19 Ag Test Letter of Authorization<sup>1</sup>, along with the authorized Fact Sheet for Healthcare Providers, the authorized Fact Sheet for Patients, and authorized labeling will be available on the FDA website post authorization: <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas>.

However, to assist clinical laboratories in using the SPERA™ COVID-19 Ag Test, the relevant Conditions of Intended Authorization are listed below:

- Authorized laboratories using your product must include with the test result reports, all Fact Sheets. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- Authorized laboratories using your product must use the product as outlined in the "authorized labeling". Deviations from the authorized procedures, including the authorized instruments, authorized clinical specimen types, authorized control materials, authorized other ancillary reagents and authorized materials required to use the product are not permitted.
- Authorized laboratories that receive your product must notify the relevant public health authorities of their intent to run the product prior to initiating testing.
- Authorized laboratories using your product must have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- Authorized laboratories must collect information on the performance of your product and report to DMD/OHT7-OIR/OPEQ/CDRH (via email: CDRH-EUA-Reporting@fda.hhs.gov) and SPERA™ Customer Support (support@xtravahealth.com) any suspected occurrence of false positive or false negative results and significant deviations from the established performance characteristics of the product of which they become aware.
- All operators using your product must be appropriately trained in performing and interpreting the results of your product, use appropriate personal protective equipment when handling this kit, and use your product in accordance with the authorized labeling.
- SPERA™ COVID-19 Ag, authorized distributors, and authorized laboratories using your product will ensure that any records associated with this EUA are maintained until otherwise notified by FDA. Such records will be made available to FDA for inspection upon request.

## Clinical Performance

The clinical performance characteristics of the SPERA™ COVID-19 Ag Test were evaluated in a multi-site prospective study in the U.S. in which patients were sequentially enrolled and tested between June 4, and July 22, 2021.

The performance of the SPERA™ COVID-19 Ag Test was established with nasal swabs collected from individual symptomatic patients (within 5 days of onset) who were suspected of COVID-19. The study was conducted at three (3) point of care (POC) sites in U.S. by 14 test operators. No training on the use of the test was provided to the operators. Operators used the Quick Reference Guide to perform testing. There were 337 subjects of the 342 total enrolled subjects (98.5%) that were determined to be eligible per the protocol (for the primary objective). A total of 337 eligible subjects were enrolled of which 49 (14.5%) were COVID positive and 288 (85.4%) were COVID negative by the comparator test. Test results were compared to the results from a highly sensitive EUA approved COVID-19 RT-PCR test.

External control testing, using the SPERA™ COVID-19 Ag Test Positive and Negative Controls, was performed prior to sample testing each day, at all study sites.

The agreement between the SPERA™ COVID-19 Ag Test and the RT-PCR comparator method are presented in the table below.

**Table 1. SPERA™ COVID-19 Ag Test Performance Against Comparator Method**

SPERA™ COVID-19 Ag Test	Comparator #1 Method		
	Positive	Negative	Total
Positive	45	9 <sup>2</sup>	54
Negative	4	279	283
Total	49	288	337
Positive Percent Agreement (PPA)	(45/49) x 100% = 91.8% (95% CI: 80.4 to 97.7%)		
Negative Percent Agreement: (NPA)	(279/288) x 100 = 96.9% (95% CI: 94.2 to 98.6%)		

## Patient Demographics

Patient demographics (gender, age, time elapsed since onset of symptoms) are available for the 337 samples used in the analysis. The table below shows the positive results stratified by patient age for the SPERA™ COVID-19 Ag Test.

**Table 2. SPERA™ COVID-19 Ag Test Positive Results by Age Group**

Age	Comparator Positive	Comparator Negative	Sum	% Positivity Rate
5 to 21 years	6	29	35	17.1
22 to 59 years	30	195	225	13.3
≥ 60	13	64	77	16.9
Sum	49	288	337	14.5

**Table 3. Positive Results Stratified by Days Post-Symptom onset**

Days Post Symptom Onset	# Specimens Tested	# Positive Specimens	Comparator Positive
0	29	2	2
1	100	8	8
2	86	10	9
3	71	17	15
4	41	15	13
5	10	2	2
Total	337	54	49

## Analytical Performance

### Limit of Detection (LOD) – Analytical Sensitivity

The limit of detection for the SPERA™ COVID-19 Ag Test was determined by using radiation-inactivated SARS-CoV-2 (isolate USA-WA1/2020) spiked onto sterile nasal swabs. A preliminary LoD was determined by first testing serial two-fold dilutions of gamma irradiated SARS-CoV-2 (isolate USA-WA1/2020) stock diluted in pooled negative nasal matrix (PNM) in triplicate which was further confirmed by an additional 20 replicates. The confirmed LoD was determined to be 1.56 x 10<sup>3</sup> TCID<sub>50</sub>/mL.

**Table 4. LoD of the SPERA™ COVID-19 Ag Test**

Virus Concentration (TCID <sub>50</sub> /mL)	Number of Positive/Total	% Detected
1.56 x 10 <sup>3</sup>	20/20	100%

### Cross-Reactivity (Analytical specificity) and Microbial Interference

Analytical specificity of the SPERA™ COVID-19 Ag Test was evaluated with a panel of sixteen (16) viruses, ten (10) bacteria, three (3) fungi, and pooled nasal wash. Final target organism concentrations were tested at ≥ 1.43 x 10<sup>5</sup> TCID<sub>50</sub>/mL, 1.0 x 10<sup>5</sup> PFU/mL, or 1.43 x 10<sup>5</sup> CEID<sub>50</sub>/mL for viruses, and ≥ 1.0 x 10<sup>6</sup> cfu/mL for bacteria and fungi.

The microbial interference was performed with the same panel of microorganisms at the same concentrations in the samples that were spiked with SARS-CoV-2 at 3X LoD. The samples were tested in triplicates for both cross-reactivity and interference studies. No cross-reactivity and no microbial interference were observed. The results for cross-reactivity and microbial interference are presented in the table below.

1 - The letter of authorization refers to, "Laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. §263a, that meet the requirements to perform high, moderate, or waived complexity tests. This test is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation." as "authorized laboratories."

2 - Following discordant analysis two false negative samples were confirmed as positive samples.

**Table 5. Cross- Reactivity and Microbial Interference Testing of the SPERA™ COVID-19 Ag Test**

Pathogen	Cross Reactivity Results	Microbial Interference Results
Human coronavirus 229E	No Cross- Reactivity	No Interference
Human coronavirus OC43	No Cross- Reactivity	No Interference
Human coronavirus NL63	No Cross- Reactivity	No Interference
SARS-coronavirus	Cross- Reactivity	Interference
MERS-coronavirus	No Cross- Reactivity	No Interference
Adenovirus	No Cross- Reactivity	No Interference
Human metapneumovirus 4 Type B2	No Cross- Reactivity	No Interference
Parainfluenza virus 1	No Cross- Reactivity	No Interference
Parainfluenza virus 2	No Cross- Reactivity	No Interference
Parainfluenza virus 3	No Cross- Reactivity	No Interference
Parainfluenza virus 4b	No Cross- Reactivity	No Interference
Influenza A	No Cross- Reactivity	No Interference
Influenza B	No Cross- Reactivity	No Interference
Enterovirus 68	No Cross- Reactivity	No Interference
Respiratory syncytial virus	No Cross- Reactivity	No Interference
Rhinovirus	No Cross- Reactivity	No Interference
Haemophilus influenzae	No Cross- Reactivity	No Interference
Streptococcus pneumonia	No Cross- Reactivity	No Interference
Streptococcus pyogenes	No Cross- Reactivity	No Interference
Candida albicans	No Cross- Reactivity	No Interference
Bordetella pertussis	No Cross- Reactivity	No Interference
Mycoplasma pneumonia	No Cross- Reactivity	No Interference
Chlamydia pneumoniae	No Cross- Reactivity	No Interference
Legionella pneumophila	No Cross- Reactivity	No Interference
Mycobacterium tuberculosis	No Cross- Reactivity	No Interference
Pneumocystis carinii	No Cross- Reactivity	No Interference
P. jirovecii-S. cerevisiae	No Cross- Reactivity	No Interference
Staphylococcus aureus subsp. aureus	No Cross- Reactivity	No Interference
Staphylococcus epidermidis	No Cross- Reactivity	No Interference
Pooled Negative Matrix	Negative	Positive

### High Dose Hook Effect

A hook effect study was completed with inactivated SARS-CoV-2 by spiking the highest concentration possible (concentration  $2.8 \times 10^6$  TCID<sub>50</sub> / ml) of SARS-CoV-2 onto swabs in triplicate and processing the replicates according to the Instructions for Use. No high dose Hook Effect was observed.

**Table 6. Hook Effect Test Results**

Test Concentration TCID <sub>50</sub> /mL	Replicates	Positive Results
$2.80 \times 10^6$	3	3/3

### Endogenous Interfering Substances

The SPERA™ COVID-19 Ag Test was evaluated for performance in the presence of potentially interfering substances that might be present in a respiratory specimen. Negative specimens were evaluated in triplicate to confirm that the substances were not cross-reactive with the test. Specimens containing SARS-CoV-2 at a concentration near the limit of detection were also evaluated in the presence of the substances in triplicate to confirm that SARS-CoV-2 could still be detected. There was no interference observed for any of the tested substances (Table 7).

**Table 7. Interfering Substances Testing of the SPERA™ COVID-19 Ag Test**

Interfering Substance	Concentration	Number Pos / # Tested with unspiked	Number Pos / # Tested with virus spiked
Human Whole Blood (EDTA tube)	4% v/v	(0/3)	(3/3)
Mucin (porcine stomach, type II)	0.5%	(0/3)	(3/3)
Chloraseptic (Menthol/Benzocaine)	1.5 mg/mL	(0/3)	(3/3)
Naso GEL (NeilMed)	5% v/v	(0/3)	(3/3)
Nasal Drops (Phenylephrine)	15% v/v	(0/3)	(3/3)
Nasal Spray (Oxymetazoline)	15% v/v	(0/3)	(3/3)
Nasal Spray (Cromolyn)	15% v/v	(0/3)	(3/3)
Zicam	5% v/v	(0/3)	(3/3)
Homeopathic (Alkalol)	10% v/v	(0/3)	(3/3)
Sore Throat Phenol Spray	15% v/v	(0/3)	(3/3)
Tobramycin	4 pg/mL	(0/3)	(3/3)
Mupirocin	10 mg/mL	(0/3)	(3/3)
Tamiflu (Oseltamivir Phosphate)	5 mg/mL	(0/3)	(3/3)
Fluticasone Propionate	5% v/v	(0/3)	(3/3)

### Symbol Glossary

	<b>In vitro diagnostic medical device</b>	Indicates a medical device that is intended to be used as an in vitro diagnostic medical device
	<b>Manufacturer</b>	Indicates the medical device manufacturer
	<b>Use by</b>	Indicates the date after which the medical device is not to be used
	<b>Batch code</b>	Indicates the manufacturer's batch code to identify the batch or lot
	<b>Catalog number</b>	Indicates the manufacturer's catalogue number to identify the medical device
	<b>Temperature limit</b>	Indicates the temperature limits to which the medical device can be safely exposed
	<b>Do not reuse</b>	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure
	<b>Caution</b>	Indicates the need for the user to consult the instructions for use for important cautionary information
	<b>Prescription use only</b>	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner
	<b>Contains sodium azide</b>	Contains sodium azide, toxic hazard addressed in the SDS document
	<b>Consult instructions for use</b>	Indicates the need for the user to consult the instructions for use
	<b>Keep Away from Sunlight</b>	Indicates a medical device that needs protection from light sources.
	<b>Do Not Use if Package is Damaged</b>	Indicates a medical device that should not be used if the package has been damaged or opened.
	<b>Sterilized Using Irradiation</b>	Indicates a medical device that has been sterilized using irradiation.






## Specimen Collection and Handling

### Preparation for Sample Collection

1. Determine the number of samples to be collected.
2. Prepare one sterile swab, one sample collection tube, one sample collection cap, and one test device for each patient sample to be collected.  
**NOTE:** If you notice the swab is damaged prior to use, discard the swab.
3. Use components immediately after opening.

### Sample Collection

Collect the sample using the provided swab following the CDC guidelines <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>






				
<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
<b>OPEN</b>	<b>INSERT SWAB</b>	<b>COLLECT SAMPLE</b>	<b>REMOVE SWAB</b>	<b>REPEAT</b>
Remove the swab from the container, being careful not to touch the soft end with your hand.	Insert the entire collection tip of the swab provided (usually 1/2 to 3/4 of an inch, or 1 to 1.5 cm) inside the nostril.	Firmly sample the nasal wall by rotating the swab in a circular path against the nasal wall at least 4 times.  Take approximately 15 seconds to collect the specimen. Be sure to collect any nasal drainage that may be present on the swab.	Gently remove the swab.	Repeat in the other nostril using the same swab.

## Specimen Transport and Storage

Specimens collected for the SPERA™ COVID-19 Ag Test should be tested as soon as possible and should not be stored or transported. The sample is stable for up to 60 minutes after preparation. Testing collected specimen after 60 minutes could yield false results.

## Test Procedure


⚠ Contains Sodium Azide ⚠ Do not use if seal is broken

				
<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>
<b>DIP SWAB</b>	<b>MIX WELL</b>	<b>REMOVE SWAB</b>	<b>ATTACH CAP</b>	<b>ADD SAMPLE</b>
Remove the aluminum seal from the sample collection tube. Immerse the swab containing the collected sample into the liquid in the sample collection tube.	Pinch the tip of the swab from the outside of the tube. Continue to pinch the outside of the tube to release the sample from the swab by moving the swab up and down while rotating. Continue the process for 15 seconds.	Pinch the outside of the tube above the tip at the top of the tube and pull out the swab (while squeezing the tube) to remove any remaining sample and suspension solution from the swab.	Securely attach the sample collection cap to the collection tube containing the sample.	Slowly turn it vertically upside down, pinch the tube, and add 3 drops to the sample well on the test device. The first drop may contain bubbles, but this will not affect the test results.  <b>Note:</b> Fewer than three drops may result in false negative results.

## IMPORTANT


- Once the sample has been applied to the sample well, allow the test device to sit undisturbed stand at 15-30°C for 15 minutes on a level surface, then interpret results. Do not read results if more than 30 minutes after applying the sample to the test device have elapsed.
- Note: Interpretation of the test before 15 minutes or after 30 minutes could yield false results. Do not interpret the test before 15 minutes or past 30 minutes.

## Interpretation of Results



Read Test Device

### POSITIVE RESULT




If a red "T" test line and a blue "C" control line are visible at or before 15 minutes, the test is positive and the result is valid.

**Note: Additional confirmatory testing with a molecular test for positive results may also be necessary, if there is a low likelihood of SARS-CoV-2 infection, such as in individuals without known exposures to SARS-CoV-2 or residing in communities with low prevalence of infection.**

#### LOOK CLOSELY!

The "T" line can be very faint. Any red/pink line visible here indicates a positive result.

### NEGATIVE RESULT



If only the blue control line under the letter "C" and no red test line under the letter "T" appears by the end of the 15-minute time frame, the sample is considered negative and results are valid.

**Note: Negative results should be treated as presumptive and may be confirmed with a molecular assay, if necessary, for patient management.**

### INVALID RESULT



If a blue control line does not appear under the letter "C" by the end of the 15-minute time frame, regardless of the status of the red test line under the letter "T", the test is invalid. Discard test materials and repeat the test.

## ORDERING AND CONTACT INFORMATION

### Technical Support

Contact us at 1-888-987-2821 or email us at [support@xtravahealth.com](mailto:support@xtravahealth.com)



Manufactured for:  
**Xtrava Inc., DBA Xtrava Health**  
3080 Olcott St. C201, Santa Clara CA 95054  
1-888-987-2821  
[support@xtravahealth.com](mailto:support@xtravahealth.com)  
[www.xtravahealth.com](http://www.xtravahealth.com)